

## Point-of-Care, Laptop/Phlebotomy Select Cart



### Features

- Tough, powder-coated, steel exterior, center cabinet
- ClintonClean™ solid HDPE, plastic base and top
- Easy-clean, plastic, 19" W x 14" D x 7.75" H elevated surface for laptop
- Top has integrated easy-grip handle and 3/4" catch channel to help prevent objects from rolling off the surface, a 4" x 5" integrated cutout to hold items and three special cutouts to hold other essential supplies
- ClintonClean™ door with clear PETG panels makes identifying supplies quick and easy
- Easy-clean, white plastic interior sides
- Three 3.5" medical, drawer bins are standard. Choose options for other medical, drawer bin configurations. All medical, drawer bins come with dividers.
- All drawer bins are on rust-free plastic tracks
- Standard 4", hospital style casters, two locking
- 15.75" x 4" glove box holder with storage compartment; two, 15.75" x 3" open, side bins; and one 15.75" x 4", clear, covered, tilt-lid, side bin
- 10.75" x 7.5" x 12" deep waste receptacle with 3/4" ClintonClean™ holder
- Keyed lock secures\* cabinet
- All parts are field replaceable, if ever damaged
- 3-year, limited warranty

**NOTE: Laptop and supplies pictured not included.**

**\* Lock not intended to secure narcotics or other restricted items.**

### SPECIFICATIONS

|                            |                       |
|----------------------------|-----------------------|
| Top                        | 32" L x 19" W         |
| Base                       | 22.5" L x 19" W"      |
| Height (less laptop shelf) | 37.5"                 |
| Cabinet                    | 17" L x 14" W         |
| 3.5" Blue Bin              | 15" L x 9" W x 3.5" D |
| 5" Blue Bin                | 15" L x 9" W x 5" D   |

### Prop 65 Warning



**WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1



**MODEL 67-145**  
Three 3.5" Bins



**Side**  
Clear, covered, tilt-lid, side bin and trash receptacle and holder



**Open Door**



**67-145 with option 304**  
Three 3.5" Bins & One 5" Bin



**67-145 with option 305**  
Four 5" Bins & One 3.5" Bin



**67-145 with option 306**  
Six 3.5" Bins

## Additional Point-of-Care, Phlebotomy Select Cart Options



-601 Electronic Lock

**NOTE:** Electronic Lock not intended to secure narcotics or other restricted items.



-602 Five inch Casters (upgrade from 4 inch)



-611 No 4"x 5" Cutout—No Charge

## Dividers in all Bins (make up to 6 sections)



### DIVIDERS INCLUDED WITH EACH BIN:

- Two Short & One Long

### ADDITIONAL DIVIDERS

---

#### CHOOSE OPTIONS:

-735 (3.5" DEEP SET)  
Three Short & One Long

-750 (5" DEEP SET)  
Three Short & One Long