

Point-of-Care, Slimline Cart



Features

- Tough, powder-coated, steel exterior, center cabinet
- ClintonClean™ solid HDPE, plastic base and top
- Top has integrated easy-grip handle and 3/4" catch channel to help prevent objects from rolling off the surface
- Tray stops on top helps secure tray when moving
- 17" width allows for maneuvering in tight spaces
- ClintonClean™ door with clear PETG panels makes identifying supplies quick and easy
- Easy-clean, white plastic interior sides
- Three 3.5" medical, drawer bins are standard. Choose options for other medical, drawer bin configurations. All medical, drawer bins come with dividers.
- All drawer bins are on rust-free plastic tracks
- Standard 4", hospital style casters, two locking
- Glove box holder and three 11.25" x 3" bins located under the handle
- Keyed lock secures* cabinet
- All parts are field replaceable, if ever damaged
- 3-year, limited warranty

NOTE: Tray and supplies pictured not included.

* Lock not intended to secure narcotics or other restricted items.

SPECIFICATIONS

Top	25" L x 17" W
Base	21.25" L x 17" W"
Height	37.5"
Cabinet	17" L x 14" W
3.5" Blue Bin	15" L x 9" W x 3.5" D
5" Blue Bin	15" L x 9" W x 5" D

Prop 65 Warning



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



MODEL 67-364
Three 3.5" Bins



Back
Glove box holder
& 3 sliding side bins



Open Door



with option 304
Three 3.5" Bins & One 5" Bin



with option 305
Four 5" Bins & One 3.5" Bin



with option 306
Six 3.5" Bins

Additional Point-of-Care, Slimline Cart Options



-601 Electronic Lock

NOTE: Electronic Lock not intended to secure narcotics or other restricted items.



-602 Five inch Casters (upgrade from 4 inch)

Dividers in all Bins (make up to 6 sections)



DIVIDERS INCLUDED WITH EACH BIN:

- Two Short & One Long

ADDITIONAL DIVIDERS

CHOOSE OPTIONS:

-735 (3.5" DEEP SET)
Three Short & One Long

-750 (5" DEEP SET)
Three Short & One Long