

# CANDo®

## Treatment Table with Adjustable Back & Shelf



### Key Features

- Available with our Quick Ship program, meaning it ships to you and your clients faster
- 2-section hardwood treatment table with an adjustable back and hardwood shelf
- 2" black upholstered foam mat top with seamless rounded corners
- Adjustable back features 5 levels of adjustment
- Measures: 78" L x 30" W; maximum capacity: 400 lbs

15-4225 CanDo® Treatment Table, with adjustable back and shelf, 78" x 30", Quick Ship