



C-Matt™ Replacement VersaCare® Foam Mattress

Top of the Line Mattress for a Fraction of the Cost

Our C-Matt Replacement VersaCare mattress is a quality patient support surface that offers exceptional safety and comfort benefits. This mattress is custom designed to meet your every need — **patient comfort and safety, infection control, and easy cleaning and maintenance.** And perhaps the best part of all, we offer all of this at a fraction of the cost of other manufacturers.



Features & Benefits

C-Matt™ Replacement VersaCare Foam Mattress technology:

- **Adjustable length.** Adjusts between 80 – 84 inches
- **Magnet technology.** Secures mattress to bedframe
- **Custom stay-put design.** Conforms to bedframe
- **Replacement covers & refresher components.** Saves money — Refresh, don't replace

Cover features include:

- **Patented seal.** Prevents fluids from penetrating mattress
- **Bi-directional stretch.** Adds to comfort and durability
- **Microclimate technology.** Absorbs and disperses body heat and perspiration
- **Full perimeter zipper.** Offers quick, easy removal
- **Antimicrobial and hypoallergenic.** Inhibits bacterial growth and allergens
- **Medical-grade tri-laminate vinyl bottom.** Extends mattress life



Foam features include:

- **Multi-zone foam core.** Enhances patient comfort and safety
- **Viscoelastic memory foam topper.** Redistributes pressure at high-risk body points
- **Oversized firm perimeter side rails.** Provides security and support
- **Shear liner.** Moisture vapor permeable and eliminates friction between foam core and cover
- **Fire barrier.** Meets all state and federal fire safety standards and Cal 129/Fed 1633



Specifications:

Specifications

Brand Name
Depth
Facility Type
Fire Barrier
Latex Free
Length
Mattress Core Material

Mattress Cover Material

Risk Level
Warranty
Weight Capacity
Width

Values

C-Matt Versa Care Foam
6"
Acute Care
Yes
Yes
80" to 84" Expandable
High Density & Visco Foam
Deep Cut Swirl Convolute
Expandable Visco Heel Section
Bi-Directional Stretch Urethane Top
Tri-Laminate Vinyl Bottom
High-Risk Early Stage 3
5 Years
500 lbs.
36"

Packaging Details

Unit of Measure
Each

Weight
37 lbs.

Package Dimensions
Palletized 10 X 40 X 90
Up to 10 Units per Pallet

Call for a quote on Casco's C-Matt Replacement VersaCare mattress and start saving today!