



C-Matt™ TotalCare Sport LAL

16 Hose Complete Mattress Cover Assembly

The fully integrated TotalCare SpO2RT Bed System Replacement Cover synchronizes with the frame when any function is activated.

Key features of this cover are:

- 1) Quality and functionality is the same as the Hill-Rom TotalCare SpO2RT Bed System Cover
- 2) Pulmonary fuchsia color
- 3) Bottom Cover completely sealed to prevent fluid ingress
- 4) Modified larger X-Ray pocket for easy access and larger photo capabilities
- 5) Full perimeter zipper for ease of service
- 6) Equipped to perform all functions as the original Hill-Rom TotalCare SpO2RT Bed System Cover
- 7) Casco Five Year Limited Warranty



Cover features include:

- 1. Durable BravoSilver (Silver Bold) Infused Polyurethane.**
 - Controlled release of silver ions resulting in a >99.9% reduction of bacteria
 - Highly resistant to stains and harsh disinfectants
 - Wipe clean between patients
- 2. Managing microclimate through our support surface cover will minimize heat build-up and reduce the likelihood of perspiration and moisture accumulation.**
 - The replacement surface provides airflow close to the patient where in is most effective
 - Provides lower surface temperatures
 - The coverlet is vented to allow moisture vapor to be transported outside the surface and away from the patient.

Low Air Loss Capabilities

Pressure redistribution

- The replacement cover fits over the multiple zones that allow the patient to immerse and the mattress uniformly envelop to patient

LAL Comfort Topper

- Provides airflow directly under the patient to keep them cool and dry while providing unsurpassed comfort to the patient.

Low Friction, low shear coverlet

- The breathable and very comfortable low friction, low shear coverlet provides superior skin protection

