

AGP34RFLCALADA



32.75" x 19" x 19.75" (H x W x D)

General purpose ADA height 2-door refrigerator-freezer with digital refrigerator control, internal fan, combination lock, and NIST calibrated alarm/thermometer

Highlights:

Combination refrigerator-freezer for general purpose storage in facilities

Digital controls allow you easily manage the refrigerator temperature between 2 and 8°C

NIST calibrated thermometer included to help monitor the interior temperature

Product Features:

Refrigerator-freezer combination	2 to 8°C refrigerator section for basic storage applications, with an upper freezer compartment ideal for ice packs
Combination lock	Factory-installed front lock secures both compartments with a user-set code. Two keys are included to override the combination as needed
Digital thermostat (refrigerator)	Easy operation with a factory-installed digital thermostat that allows you to manage the refrigerator temperature between 2 and 8°F (NOTE: the freezer temperature will self-adjust based on the refrigerator setting)
Alarm/thermometer included	Temperature monitoring device (TMD) displays the current and min/max temperature of the refrigerator interior, with an audible alarm to alert if the unit goes out of range
Buffered temperature probe	Thermometer probe is encased in a glycol-filled bottle to better simulate the temperature of stored contents instead of refrigerator air temperature
NIST calibrated external temperature readout	Thermometer is NIST calibrated in our ISO/IEC 17025:2017 accredited laboratory and ships with a 2-year certificate of calibration for your records
Fan-forced cooling	Refrigerator section includes a factory-installed fan to help ensure a more evenly cooled interior
Automatic defrost refrigerator	Lower refrigerated section does not require manual defrosting

AGP34RFLCALADA Specifications:

Overview	
Nominal Height	32.0" (81 cm)
Height of Cabinet	32.63" (83 cm)
Height to Hinge Cap	32.75" (83 cm)
Nominal Width	19.0" (48 cm)
Width	19.0" (48 cm)
Nominal Depth	20.0" (51 cm)
Depth	19.75" (50 cm)
Depth with door at 90°	36.75" (93 cm)
Capacity	3.2 cu.ft. (91 L)
Defrost Type	Cycle
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	1.1
Voltage/Frequency	115 V AC/60 Hz
Weight	55.0 lbs. (25 kg)
Parts & Labor Warranty	1 Year
Compressor Warranty	5 Years
UPC	761101086033
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Yes
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	18.5" (47 cm)
Refrigerator - Interior Width	15.75" (40 cm)
Refrigerator - Interior Depth	15.0" (38 cm)
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	2
Refrigerator - Full Door Shelves	2
Refrigerator - Temperature Range	2 to 8°C
Freezer - Interior Height	7.5" (19 cm)
Freezer - Interior Width	14.0" (36 cm)
Freezer - Interior Depth	14.0" (36 cm)

Freezer - Full Door Shelves	1
Thermostat Type	Digital
Fan Type	Interior
Refrigerant Type	R600a
Refrigerant Amount	1.16 oz.
High Side PSI	250.0
Low Side PSI	88.0
Level Legs	2
Compressor Step Height	9.0" (23 cm)
Compressor Step Width	15.75" (40 cm)
Compressor Step Depth	5.0" (13 cm)

