

# ARFG15PV



74.63" x 23.5" x 25" (H x W x D)

PHARMA-VAC Performance Series combination 2-door/2-compressor refrigerator-freezer for pharmacy, medical, and life science applications

### Highlights:

Purpose-built design to store medication, pharmaceuticals, and other sensitive materials

Independent dual compressor refrigerator and freezer storage in one slim-fitting footprint to maximize storage capacity without compromising performance

Intelligent microprocessor digital temperature controller for optimum temperature management

### Product Features:

Performance Series Refrigerator-Freezer Combination	Purpose-built all-refrigerator and all-freezer for vaccine storage in medical, pharmaceutical, and life science applications
Purpose-built design	Designed and purpose-built for pharmacy, medication, and vaccine applications to support meeting CDC/VFC/VFA guidelines
+/-1°C Variation	The published temperature variation is derived from the maximum deviation of an NTC sensor in a 1 oz. vial located nearest the chamber geometric center during a 24 hour test period
Refrigerator window	Triple paned glass on the upper refrigerator section provides a clear view of interior contents
Space-saving design	Allows you to create a full refrigerator-freezer with independent controls in a slim-fitting footprint
Antimicrobial silver ion-coated handles	Powder-coated handles with naturally occurring silver ions to help reduce the spread of germs
Adjustable temperature	Intuitive controls let you adjust the interior temperature between +2 to +8°C (refrigerator) and -30° to -10°C (freezer)
Temperature alarms	Visual and audible alarms alert if the interior temperature rises or falls out of set range
Remote alarm contacts	Electronic controller is equipped with a volt-free contact for alarm forwarding to an external

## ARFG15PV Specifications:

Overview	
Height of Cabinet	74.63" (190 cm)
Width	23.5" (60 cm)
Width with Door Open	24.0" (61 cm)
Depth	25.0" (64 cm)
Depth with Handle	26.5" (67 cm)
Depth with door at 90°	46.75" (119 cm)
Capacity	10.74 cu.ft. (304 L)
Door	White
Cabinet	White
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	2.0
Voltage/Frequency	115 V AC/60 Hz
Weight	205.0 lbs. (93 kg)
Parts & Labor Warranty	2 Years
Compressor Warranty	5 Years
UPC	761101156705
Refrigerator-Freezer	
Door Swing	RHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Interior Light	Yes
Refrigerator - Interior Height	33.0" (84 cm)
Refrigerator - Interior Width	19.5" (50 cm)
Refrigerator - Interior Depth	17.5" (44 cm)
Refrigerator - Capacity	6.74 cu.ft
Refrigerator - Shelf Type	Wire
Refrigerator - Shelf Qty	3
Refrigerator - Temperature Range	2 to 8°C
Refrigerator - Refrigerant Type	R600a
Refrigerator - Refrigerant Amount	1.3 oz.
Freezer - Interior Height	17.75" (45 cm)
Freezer - Interior Width	15.25" (39 cm)
Freezer - Interior Depth	16.5" (42 cm)

Freezer - Capacity	4.0 cu.ft
Freezer - Shelf Type	Wire
Freezer - Shelf Qty	2
Thermostat Type	Digital
Fan Type	Interior
Wheels	4
Freezer - Temperature Range	-30 to -10°C
Freezer - Refrigerant Type	R290a
Freezer - Refrigerant Amount	1.2 oz.
Sound Level	44.5 dB

