


# PARALLEL BARS

## PLATFORM MOUNT

### STANDARD FEATURES:

- Load capacity: 500 lbs.
- Bars: 19-30"W and 26-44"H
- Platform 48"W
- 19-30"W x 84" or 120"L x 26-44"H
- Stainless steel hand rails (1.5" Diameter)
- Posi-Lok locking knobs
- Some assembly required



AM-710

### STANDARD CONFIGURATIONS:

**AM-707** - 7' Platform Mount Parallel Bars  
19-30"W x 84"L x 26-44"H - Platform 48"W

**AM-710** - 10' Platform Mount Parallel Bars  
19-30"W x 120"L x 26-44"H - Platform 48"W  
2 Piece Platform

### OPTIONS:

**AM-714** - Abduction Board for Parallel Bars  
.75"W x 94"L x 5"H

**AM-715** - Balance Beam for Parallel Bars  
5"W x 80"L x 2.5"H

**AM-716** - Foot Placement Ladder for Parallel Bars  
24"W x 84"L x 4"H