


# PARALLEL BARS

## FLOOR MOUNT

### STANDARD FEATURES:

- Load capacity: 500 lbs
- 84" or 120"L x 26-44"H
- Stainless steel hand rails (1.5" Diameter)
- Posi-Lok locking knobs
- Mounting hardware not included
- Some assembly required



AM-711

### STANDARD CONFIGURATIONS:

**AM-711** - 7' Floor Mount Parallel Bars  
84"L x 26-44"H

**AM-712** - 10' Floor Mount Parallel Bars  
19-30"W x 120"L x 26-44"H

### OPTIONS:

**AM-716** - Foot Placement Ladder for Parallel Bars  
24"W x 84"L x 4"H