

# PARALLEL BARS

## ELECTRIC PLATFORM MOUNT

### STANDARD FEATURES:

- Load capacity: 500 lbs
- Bars: 20-34"W and 26-54"H
- Platform 48"W
- Stainless steel hand rails (1.5" Diameter)
- Digital display on conveniently located, user-friendly control panel shows exact bar height
- Memory buttons allow storage and recall of four separate bar height settings
- Microprocessor control unit ensures synchronous operation of all four columns
- Height monitoring functions ensure reliable operation and safe switch-off in case of a fall
- Slip-resistant platform
- ETL Approved
- Some assembly required



AM-TM702 with End Control

### STANDARD CONFIGURATIONS:

<p><b>AM-TM702</b> - 7' Electric Platform Mount Parallel Bars (Single Platform) - Platform: 48"W 20-34"W x 84"L x 26-54"H</p>	<p><b>AM-TM700</b> - 10' Electric Platform Mount Parallel Bars (2 Piece Platform) - Platform: 48"W 20-34"W x 120"L x 26-54"H</p>	<p><b>AM-TM703</b> - 15' 2PC Electric Platform Mount Parallel Bars (3 Piece Platform) - Platform: 48"W 20-34"W x 180"L x 26-54"H</p>
-----------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------

### OPTIONS:


**AM-714** - Abduction Board for Parallel Bars  
.75"W x 94"L x 5"H

**AM-715** - Balance Beam for Parallel Bars  
5"W x 80"L x 2.5"H

**AM-716** - Foot Placement Ladder for Parallel Bars  
24"W x 84"L x 4"H



Close Up of Section Attachments



View of Base

