

# PARALLEL BARS

## ELECTRIC PLATFORM MOUNT

### STANDARD FEATURES:

- Load capacity: 800 lbs (400 lbs per bar)
- 8', 10', and 12' Lengths available
- Platform: 50"W (19.125"-30.5"W between columns)
- Heavy duty, rugged, patented height and width adjustment system
- Steel, multi-section platform with gray anti-bacterial and anti-microbial powder coat base
- Brushed stainless steel hand rails (2" Diameter)
- Battery powered
- Microprocessor control unit ensures synchronous operation of all four stanchions
- Height monitoring functions ensure reliable operation and safe switch-off in case of a fall
- Height and width can move independently or in unison
- Slip-resistant platform
- ETL Approved
- Some assembly required



### STANDARD CONFIGURATIONS:

**EPB10080** - 8' Electric Platform Mount Parallel Bars (3 Piece Platform)  
19.125-30.5"W x 96"L x 26.5-42.625"H

**EPB10010** - 10' Electric Platform Mount Parallel Bars (4 Piece Platform)  
19.125-30.5"W x 120"L x 26.5-42.625"H

**EPB10012** - 12' Electric Platform Mount Parallel Bars (4 Piece Platform)  
19.125-30.5"W x 144"L x 26.5-42.625"H

EPB10010