

WEIGHT SLED


STANDARD FEATURES:

- Load capacity: 250lbs
- 20"W X 36"L (Base size)
- T-Handle adjusts from 35-41"
- Hardwood construction
- Removable weight holder
- Durable plastic runners on the underside of the sled reduce floor wear and provide smooth glide performance, making the unit ideal for use in therapy, rehabilitation, and recovery environments
- Push, pull, or drag this sled to help condition and strengthen the body

SIZES AVAILABLE:

WCWS

20"W x 28"L x 33"-45"H



WCWS