

DELUXE MOBILE LEAD SHIELD WITH T-BASE

PRODUCT INFORMATION



PHILLIPS SAFETY
PRODUCTS INC.

PRODUCT INFORMATION
DELUXE MOBILE LEAD SHIELD WITH T-BASE






The Deluxe Mobile Lead Shield with T-Base is an accommodating radiation shielding for both children and adults. Made of durable Nylon fabric, it is available with 0.5mm, 1.0mm, and 2.0mm PbEq protection. This Phillips Safety lead shield is excellent for patient and staff protection during imaging procedures, due to its Posi-Lock height adjustable arm with spring-loaded pin that adjusts from 8" to 60". Plus, due to its all steel T-Base, it provides excellent stability and confidence in no-tipping due to its and allows for easy positioning around equipment and wall stands. The Deluxe Mobile Lead Shield measures 24" W x 30" H.

SPECIFICATION

DIMENSIONS: 24" H x 30" W	FEATURES: Adjustable shield height from 8" to 60"
MATERIAL: Steel T-Base	PROTECTION LEVEL: 0.5mm Pb Eq, 1.0mm Pb Eq, 2.0mm Pb Eq

**Deluxe Mobile Lead Shield
with T-Base**



LB-MB-3024
SCALE: 1/2" = 1'