

MAX MOBILE LEAD SHIELD WITH H-BASE

PRODUCT INFORMATION



PRODUCT INFORMATION
MAX MOBILE LEAD SHIELD WITH H-BASE



The Max Mobile Lead Shield with H-Base is a durable lead shield with 0.5mm PbEq or 1.0mm PbEq. With all steel H-base for maximum stability, this Phillips Safety features hardened chrome steel post and a Posi-Lock height adjustable arm with spring-loaded pin. With an open base design, it has base overall dimensions of 24" W x 24" D, and base inside dimensions of 20" W x 16" D. In addition, this lead shield is available in two different sizes – 36" x 36" and 48" x 48" – and features a heavy duty 4-caster base.

SPECIFICATION

OVERALL DIMENSIONS: 20" W x 16" D

BASE OVERALL DIMENSIONS: 24" W x 24" D

AVAILABLE SIZES: 36" x 36", 48" x 48"

PROTECTION LEVEL: 0.5mm Pb Eq, 1.1mm Pb Eq

FEATURES: H-Base, heavy duty 4-caster base, hardened chrome steel post, Posi-Lock height adjustable arm with spring-loaded pin

PRODUCT INFORMATION