



C-Matt™ Gynnie® Birthing Bed Mattress

A Quality Replacement System for the Specialty Birthing Bed Stretcher



A Quality Birthing Bed Replacement Surface

Fits Gynnie® Stretcher

Breathable, Bi-Directional Stretch Fabric

Fully Sealed Cover

4in Thick Comfort Mattress

Entire Surface or Separate Sections Available

Optional Fire Barrier for Regulatory Compliance

Top Quality Replacement Surfaces from Casco

Specifically Designed to Fit Gynnie® Stretcher

- Hook and loop placement and snap hardware included

Wipeable Cover

- Fully sonically sealed to protect foam core
- Anti-microbial and hypo-allergenic
- Easy to clean, stain resistant and durable

Surface Foam Core

- High density, anti-microbial, medical grade foam
- 4in thick Comfort Mattress

Regulatory Requirements

- Meets CAL117
- Optional fire barrier meets CAL129/FED1633

Fully sealed with sonic welding to prevent fluid ingress



Hidden zipper and placement indicators for a secure fit



Quality replacement surfaces at a 30-50% cost savings

