

C-Matt™ StatiCare Replacement Stretcher Mattress

Top of the Line Stretcher Pads at Manufacturer Direct Price

C-Matt replacement stretcher pads are quality patient support surfaces that offer exceptional safety and comfort benefits. Our pads are custom designed to meet your every need – patient comfort and safety, infection control, and easy cleaning and maintenance. All C-Matt replacement pads and mattresses are made in the USA by a Certified Woman-Owned Small Business..

What the C-Matt StatiCare Replacement Stretcher Mattress Can Do For You

C-Matt StatiCare replacement mattress is a high performing static air surface. It uses interconnected air chambers that allow the mattress to automatically adjust to patient movement. This quiet, non-powered technology provides even body weight distribution. StatiCare's multi-zone convoluted foam core and many innovative technologies improve patient comfort and reduce the risk for skin breakdown.



Features & Benefits

Cover features include:

- **Patented seal.** Radio frequency (RF) welding technology and our patented “waterfall” zipper flap prevent fluids from penetrating the mattress.
- **Bi-directional stretch.** Breathable, two-way stretch fabric provides added comfort and durability. Also prevents “hammocking,” allowing for easier patient transfer.
- **Microclimate technology.** Polyurethane-coated or polycarbonate-coated fabric absorbs body heat and perspiration, and disperses it away from the body. This keeps patients cool and dry, reducing the incidence of skin breakdown.
- **Antimicrobial and hypoallergenic.** Our fabric is treated to inhibit bacterial growth and has proven to be hypoallergenic.
- **Quick, easy cleaning.** A full perimeter zipper allows easy access to inspect the core or change the cover, if necessary.
- **Non-slip bottom.** The bottom of the pad is covered in a non-slip fabric and has hook and loop positioners to keep the pad securely in place.
- **Fabric options.** Your choice of polyurethane-coated or polycarbonate-coated fabric cover.
- **Cover options.** Replacement covers are available.
- **Polycarbonate cover option.** Designed to better withstand the harsh nature of certain disinfectants that are used to clean and sterilize in many facilities.

Foam features include:

- **Interconnected air cells.** Seven anatomically located air cells form interconnected air chambers that quietly adjust to patient movement and evenly redistribute body pressure. No electricity or noisy pump required.
- **Size and shape options.** All mattresses have two tapered corners and are 5" thick and 76" long. Width options are 26" or 29."
- **High density, antimicrobial foam bladder tray.** The one-inch custom-designed topper foam adds comfort while providing superior pressure redistribution characteristics.
- **Industry-leading safety standards.** Foam meets Cal 117-2013 fire safety standards. The optional fire barrier is in compliance with all state and federal fire safety standards. Meets Cal 129 / Fed 1633 regulations.

Specifications:

Specifications

Brand name:

Thickness:

Width:

Length:

Weight capacity:

Values

C-Matt

5"

26" or 29"

76"

500 lbs.

Latex free

