

## Casco First Responder RescueAll™ Transport, Reusable

---

The Casco First Responder RescueAll™ Transport is made with an orange antimicrobial, 14oz. tri-laminate vinyl with a black polypropylene webbing along the sides, making the product not only easy to use, but easy to clean and reuse. The product design enables it to be used as a slide, serving as an alternative means of patient transport. The Casco First Responder RescueAll™ Transport comes in an orange pouch for easy storage.



- The Casco First Responder RescueAll™ Transport is rated to transport up to 1000 pounds (454 Kg).
- The Casco First Responder RescueAll™ Transport is recommended for use by professionals for lifting and transporting patients.
- Multiple handles for ease of team/crew handling.
- Unique handle design for ease of movement.
- Note: Rescue/All is not designed for use of spinal/traction stabilization.



# Casco First Responder RescueAll™ Transport, Reusable

| <u>Product Codes</u>                                 | Part #<br>GTIN<br>UPC                    |
|------------------------------------------------------|------------------------------------------|
| Casco First Responder™ RescueAll Transport, Reusable | H12720<br>00810125100383<br>810125100383 |
|                                                      |                                          |

## Usage Instructions:

- Do Not Machine Wash or Dry
- Avoid Contact with Sharp Objects
- Clean Only with a Mild Cleaner/Disinfectant
- Do Not Use if Punctured, torn, or Excessively worn
- Only Use the Handles for transporting the patient; grip securely
- Minimum 4 persons to lift and use
- Grab a minimum of four (4) handles on each side to ensure safe weight distribution, and prevent sagging.
- The storage pouch can be hung or fastened in a vehicle or industrial facility in a convenient location for easy access.

