



C-Matt™ Bariatric 500 & 900

An Economic Choice to Meet Your Facility's Health, Safety and Patient Comfort Needs



A Quality Bariatric Replacement Surface

Two Maximum Patient Weight Limit Options

500lb and 900lb

Replaceable Cover Available in Two Fabric Styles

Many Sizing Options Available

Optional Fire Barrier for Regulatory Compliance

Top Quality Replacement Surfaces from Casco

Replaceable Cover

- Available in bi-directional stretch, breathable polyurethane fabric or leak proof nylon top fabric
 - Antimicrobial and hypo-allergenic
 - Easy to clean, stain resistant and durable
 - Patented zipper-flap design to protect surface core
 - Full perimeter zipper for easy maintenance
- 5 year Warranty**

Surface Foam Core

- Multi-zoned, high density, antimicrobial, medical grade foam
- Medium-firm support

Regulatory Requirements

- Meets CAL117
- Optional fire barrier meets CAL129/FED1633


Multiple Sizing Options

- 7 x 39/42/48/54 x 80/84
- Custom sizes available

Full perimeter zipper for easy maintenance



Unique convolute to increase air flow and patient comfort



Quality replacement surfaces at a 30-50% cost savings.

