


Hemocytometer

Our Hemocytometer is specially designed for counting erythrocytes and leukocytes. Double rulings are visible through a wide range of illumination. Neubauer rulings are 0.1mm below cover glass. The H-shaped moat offers two counting areas. Fused rhodium on glass provides distinct lines.

FEATURES

- High light transmission
- Brightly illuminated line
- Two special counting chamber cover glasses
- Light transmission between 29% and 32%



Hemocytometer

Model #	Description
CTL-HEMM-GLDR	Hemocytometer, with cover glass

