

Pipette Shakers

LW Scientific pipette shakers easily achieve consistent results in specimen separation and are especially effective for blood preparation. The LWS 6 place shaker has the option of either continuous motion or timed motion using the 15 minute timer. Compact size and durability make the LW Scientific pipette shaker convenient for any laboratory.

SPECIFICATIONS

- Portable and compact
- Capacity – 6 places
- Holds 5ml - 7ml pipettes
- Fixed speed - 2500 rpm
- 15-minute timer
- Suction-cupped feet
- Space-saving design
- LED indicator light
- Durable construction
- 110v or 220v
- One year warranty



Model #	Description
SHL-PPF7-06F1	Pipette shaker 6 places

