

# Water Bath

The LW Water Bath creates a constant temperature with an open reservoir and immersion circulator, ideal for heating temperature-sensitive samples. The stainless steel reservoir is seamless, eliminating the possibility of rust and leaks, and the rounded corners are easy to clean. The control unit achieves the desired temperature quickly and maintains a constant temperature to prevent sample overheating. Select desired setting with the touch-pad digital controller, and view the actual temperature on the red LED display. The compact design and stainless steel reservoir make this product a must for every lab.



## TANK

Capacity 10 L or 20L  
Stainless steel reservoir, corrosion  
resistant

## TEMPERATURE SETTING

Digital

## TEMPERATURE READOUT

LED

## TEMPERATURE CONTROL

PID

## TEMPERATURE SENSOR

100 ohm Platinum RTD

## TEMPERATURE RANGE

Ambient + 5° to 100° C

## STABILITY

+/- 0.2

## POWER SUPPLY

110V (220V available)

## WEIGHT

10L: 20 lbs. (9.07 kb)  
20L: 21.2 lbs. (9.6 kg)

## ACCESSORIES

Lid (flat or gabled)

## Water Bath

| Model #       | Description |
|---------------|-------------|
| WBL-10LC-SSD1 | 10 Liter    |
| WBL-20LC-SSD1 | 20 Liter    |

- Excellent Temperature Control
- Convenient Large Bath Opening
- Stainless Steel Reservoir

## INTERNAL TANK DIMS

### 10L

Height: 6" (150mm)  
Width: 9.4" (238mm)  
Depth: 11.4" (290mm)

### 20L

Height: 6" (150mm)  
Width: 19.3" (490mm)  
Depth: 11.4" (290mm)

## DIMS

### 10L

Height with Controller: 13.5" (343mm)  
Height: 7" (178mm)  
Width: 11" (275mm)  
Depth: 13.5" (343mm)

### 20L

Height with Controller: 13.5" (343mm)  
Height: 7" (178mm)  
Width: 21.5" (546mm)  
Depth: 13.5" (343mm)

