


# MODEL DHK4-49GD

46.62 CUBIC FOOT LAB & PHARMACY REFRIGERATOR


|                     | WIDTH    | DEPTH    | HEIGHT   |  | WIDTH  | DEPTH  | HEIGHT |
|---------------------|----------|----------|----------|--|--------|--------|--------|
| INTERIOR DIMENSIONS | 52.75 IN | 26.25 IN | 59.5 IN  |  | 134 CM | 67 CM  | 152 CM |
| EXTERIOR DIMENSIONS | 57.25 IN | 34.75 IN | 78.75 IN |  | 145 CM | 88 CM  | 200 CM |
| CRATED DIMENSIONS   | 85 IN    | 41 IN    | 86 IN    |  | 215 CM | 104 CM | 218 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

## SIDE VIEW



## FRONT VIEW



# SO-LOW



**NET WEIGHT : 750 LBS.**  
**SHIPPING WEIGHT: 800 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )  
INSIDE DELIVERY ( WHITE GLOVE DELIVERY )