


# MODEL MV4-2UCRDA

## 2.64 CUBIC FOOT UNDERCOUNTER REFRIGERATOR


|                     | WIDTH    | DEPTH   | HEIGHT   | WIDTH | DEPTH | HEIGHT |
|---------------------|----------|---------|----------|-------|-------|--------|
| INTERIOR DIMENSIONS | 17.5 IN  | 17.5 IN | 21.75 IN | 44 CM | 44 CM | 54 CM  |
| EXTERIOR DIMENSIONS | 21.25 IN | 22 IN   | 30 IN    | 54 CM | 56 CM | 76 CM  |
| CRATED DIMENSIONS   | 32 IN    | 32 IN   | 46 IN    | 81 CM | 81 CM | 116 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

### SIDE VIEW



### FRONT VIEW



# SO-LOW



**NET WEIGHT : 125 LBS.**  
**SHIPPING WEIGHT: 150 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )  
INSIDE DELIVERY ( WHITE GLOVE DELIVERY )