


# MODEL MV4-6UCRGDDA

## 4.6 CUBIC FOOT UNDERCOUNTER REFRIGERATOR


|                     | WIDTH    | DEPTH    | HEIGHT   | WIDTH | DEPTH | HEIGHT |
|---------------------|----------|----------|----------|-------|-------|--------|
| INTERIOR DIMENSIONS | 21.5 IN  | 20 IN    | 23.25 IN | 55 CM | 51 CM | 59 CM  |
| EXTERIOR DIMENSIONS | 25.75 IN | 24.75 IN | 32 IN    | 65 CM | 63 CM | 81 CM  |
| CRATED DIMENSIONS   | 32 IN    | 32 IN    | 46 IN    | 81 CM | 81 CM | 116 CM |

NOTE THAT SHELVES ARE SHOWN IN DEFAULT CONFIGURATION, AND CAN BE ADJUSTED IN 1" INCREMENTS

### SIDE VIEW



### FRONT VIEW



# SO-LOW



**NET WEIGHT : 140 LBS.**  
**SHIPPING WEIGHT: 165 LBS.**

### SHIPPING OPTIONS

DOCK-TO-DOCK DELIVERY ( REQUIRES FORKLIFT OR PALLET JACK )  
DOCK-TO-DOCK DELIVERY WITH LIFTGATE ( REQUIRES FORKLIFT OR PALLET JACK )  
INSIDE DELIVERY ( WHITE GLOVE DELIVERY )