

## MODEL DHK4-14SD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

|                                 |                     |
|---------------------------------|---------------------|
| <b>Capacity [ Cu. Ft. ]</b>     | <b>14 Cu. Ft.</b>   |
| <b>Capacity [ Liters ]</b>      | <b>395 Liters</b>   |
| <b>Temperature Range [ °C ]</b> | <b>2°C to 8°C</b>   |
| <b>Temperature Range [ °F ]</b> | <b>35°F to 46°F</b> |

|                          |                        |
|--------------------------|------------------------|
| Dimensions               | Width x Depth x Height |
| Exterior [ Inches ]      | 25.75" x 26.5" x 78.5" |
| Interior [ Inches ]      | 23" x 21" x 53.25"     |
| Exterior [ Centimeters ] | 65 x 67 x 199          |
| Interior [ Centimeters ] | 58 x 53 x 135          |

|                 |                                |
|-----------------|--------------------------------|
| Compressor      | Hermetic Compressor            |
| Air Circulation | Forced Air                     |
| Defrost Method  | Automatic                      |
| Refrigerant     | R600a                          |
| Exterior        | Steel                          |
| Interior        | Polystyrene ( HIPS )           |
| Insulation      | US EPA and SNAP approved       |
| Lockable Door   | Yes, Keyed                     |
| Porthole        | Yes, 1" Diameter               |
| Casters         | 4, with Leveling Feet          |
| Storage         | 6 Adjustable Shelves, 1 Basket |
| Shipping Weight | 250 lbs.                       |


#### Electrical

|                         |                            |
|-------------------------|----------------------------|
| Voltage:                | 115 Volts 60 Hertz 1 Phase |
| Amperage Line   Running | 15 Amp Dedicated   3.77A   |
| Power Consumption       | 3.71 kW/d                  |
| NEMA Configuration      | NEMA 5-15                  |
| Certification(s)        | UL Listed                  |

#### Controller

|                       |                                                                                       |
|-----------------------|---------------------------------------------------------------------------------------|
| Display Type          | Microprocessor, Digital                                                               |
| Alarms                | High/Low Temp, Door Ajar<br>Power Failure, Low Battery<br>Sensor Failure, USB Failure |
| Data Logging          | Yes, Adjustable Intervals                                                             |
| Data Download         | Yes, via USB, PDF Format                                                              |
| Min / Max Temperature | Yes, Display and Reset                                                                |
| Alarm Relay           | Yes, Dry Contacts                                                                     |
| Backup Battery        | Yes, Only Powers Alarm                                                                |

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 2 Year Parts & Labor