

MODEL DH-48RF

LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

	REFRIGERATOR	FREEZER
Total Capacity [Cu. Ft.]	21.01 Cu. Ft.	21.01 Cu. Ft.
Capacity [Liters]	595 Liters	595 Liters
Temperature Range [°C]	2°C to 8°C	0°C to -23°C
Temperature Range [°F]	35°F to 46°F	32°F to -10°F

Dimensions	Width x Depth x Height	
Interior [Inches]	23" x 27" x 55"	23" x 27" x 55"
Interior [Centimeters]	58 x 68 x 139	58 x 68 x 139
Exterior [Inches]	52" x 33.75" x 81"	
Exterior [Centimeters]	132 x 85 x 205	

Compressor	Hermetic	Hermetic
	Compressor	Compressor
Air Circulation	Forced Air	Forced Air
Defrost Method	Automatic	Automatic
Porthole	1/2" Diameter	1/2" Diameter
Refrigerator Storage	5 Adjustable Shelves	
Freezer Storage	5 Adjustable Shelves	
Refrigerant	R290	
Exterior	Stainless Steel	
Interior	Stainless Steel	
Insulation	Polyurethane	
Lockable Door	Lockable with Padlock	
Casters	4, Locking	
Shipping Weight	300 lbs.	

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line	15 Amp Dedicated
Running Amperage	6.0A Refrigerator 5.7A Freezer
NEMA Configuration	NEMA 5-15
Plugs	1 Freezer, 1 Refrigerator
Certification(s)	UL Listed

Controller

Type	Microprocessor
	Two Digital Displays
Alarms	High/Low Temp
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts



Warranty 3 Year Parts & Labor.
Additional 4 Years on compressor parts.