

MODEL DH-70RF

LAB AND PHARMACY REFRIGERATOR FREEZER COMBO

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

	REFRIGERATOR	FREEZER
Total Capacity [Cu. Ft.]	19.62 Cu. Ft.	44.85 Cu. Ft.
Capacity [Liters]	555 Liters	1270 Liters
Temperature Range [°C]	2°C to 8°C	0°C to -23°C
Temperature Range [°F]	35°F to 46°F	32°F to -10°F

Dimensions	Width x Depth x Height	
Interior [Inches]	22" x 28.5" x 55"	47" x 28.5" x 55"
Interior [Centimeters]	55 x 72 x 139	119 x 72 x 139
Exterior [Inches]	75" x 33.75" x 81.5"	
Exterior [Centimeters]	190 x 85 x 207	

Compressor	Hermetic	Hermetic
	Compressor	Compressor
Air Circulation	Forced Air	Forced Air
Defrost Method	Automatic	Automatic
Porthole	1/2" Diameter	1/2" Diameter
Refrigerator Storage	5 Adjustable Shelves	
Freezer Storage	10 Adjustable Shelves	
Refrigerant	R290	
Exterior	Stainless Steel	
Interior	Stainless Steel	
Insulation	Polyurethane	
Lockable Door	Lockable with Padlock	
Casters	6, Locking	
Shipping Weight	300 lbs.	

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line	15 Amp Dedicated
Running Amperage	10.6A Refrigerator 9.4A Freezer
NEMA Configuration	NEMA 5-15
Plugs	1 Freezer, 1 Refrigerator
Certification(s)	UL Listed

Controller

Type	Microprocessor
	Two Digital Displays
Alarms	High/Low Temp
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts



Warranty 3 Year Parts & Labor.
Additional 4 Years on compressor parts.