

MODEL DHH4-17SDFMS

FLAMMABLE MATERIAL STORAGE REFRIGERATOR 1°C TO 10°C

Capacity [Cu. Ft.]	17 Cu. Ft.
Capacity [Liters]	481 Liters
Temperature Range [°C]	1°C to 10°C
Temperature Range [°F]	33°F to 50°F

Dimensions	Width x Depth x Height
Exterior [Inches]	30.75" x 28.5" x 67.75"
Interior [Inches]	26.25" x 20.75" x 58.5"
Exterior [Centimeters]	76 x 72 x 172
Interior [Centimeters]	66 x 52 x 148

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Steel
Insulation	High Density Urethane
Lockable Door	1 Door, Keyed Lock
Door Material	Insulated Steel
Porthole	3/4" Diameter
Leveling Legs	4
Storage	5 Adjustable Shelves
Shipping Weight	240 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 3A
NEMA Configuration	NEMA 5-15
Certification(s)	ETL Listed

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



» Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.

Warranty 2 Years Parts & Labor.