

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



MEDICAL RESEARCH | LABORATORY STORAGE | PHARMACY USE | NON-CRITICAL APPLICATIONS

Economic Storage


Our general purpose refrigerators maintain temperatures between 2°C to 8°C. These units are a less expensive option for general storage of non-critical items; or products which do not require an alarm or precise temperature tolerances.

Designed for Universal Applications

Our general purpose units come standard with a digital microprocessor control, door ajar alarm, power failure alarm, low battery sensor failure, and USB failure sensor. These units also come with a modern display for ease of use, and data logging with PDF download capability.

Electrical Certification

All general purpose models are UL Listed.



MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
--------------	------------	--------	-------------------	------------------------------------	------------------------------------	--------------	-------------

GLASS DOOR UNDERCOUNTER

MV4-2UCRGDDA	2.64	75	2°C to 8°C	17.5" x 17.5" x 21.75"	21.25" x 22" x 30"	AUTOMATIC	125 LBS
MV4-6UCRGDDA	4.6	130	2°C to 8°C	21.5" x 20" x 23.25"	25.75" x 24.75" x 32"	AUTOMATIC	140 LBS

SOLID DOOR UNDERCOUNTER

MV4-2UCRDA	2.64	75	2°C to 8°C	17.5" x 17.5" x 21.75"	21.25" x 22" x 30"	AUTOMATIC	125 LBS
MV4-6UCRDA	4.6	130	2°C to 8°C	21.5" x 20" x 23.25"	25.75" x 24.75" x 32"	AUTOMATIC	140 LBS

GLASS DOOR UPRIGHT

DHK4-10GD	11.12	315	2°C to 8°C	22.75" x 21" x 44"	25.75" x 26.5" x 69.5"	AUTOMATIC	215 LBS
DHK4-23GD	23	650	2°C to 8°C	23.75" x 28.5" x 59.75"	28.25" x 35" x 78.25"	AUTOMATIC	400 LBS
DHK4-49GD	46.62	1320	2°C to 8°C	52.75" x 26.25" x 59.5"	57.25" x 34.75" x 78.75"	AUTOMATIC	750 LBS

TOUCHSCREEN GLASS DOOR UPRIGHT

DHK4-23GD-T	23	650	2°C to 8°C	23.75" x 28.5" x 59.75"	28.25" x 35" x 78.25"	AUTOMATIC	400 LBS
DHK4-2349GD-T	46.62	1320	2°C to 8°C	52.75" x 26.25" x 59.5"	57.25" x 34.75" x 78.75"	AUTOMATIC	750 LBS

SOLID DOOR UPRIGHT

DHK4-14SD	14	395	2°C to 8°C	23" x 21" x 53.25"	25.75" x 26.5" x 78.5"	AUTOMATIC	215 LBS
-----------	----	-----	------------	--------------------	------------------------	-----------	---------