

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



LABORATORY RESEARCH & TESTING | PHARMACY STORAGE | GENERAL PURPOSE REFRIGERATION

Pass-Through Convenience

So-Low Pass-Through Refrigerators provide the reliability and security essential for the protection of sensitive lab samples in cleanrooms, and environmentally controlled rooms. Ideal for demanding requirements of daily lab applications, our range of refrigerators offers the dependability and safety required for your daily lab applications.

Humidity & Stability Precision

Our Stability and Humidity Chambers are manufactured with a heavy duty double wall Interior and exterior. The cabinet and door are insulated and sealed with 100% UL Class 1 rigid polyurethane foam insulation, reinforcing the precision of stability within the chamber.

Electrical Listing

Pass through units are UL Listed. Humidity & Stability chambers are ETL listed.

MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
-----------------	---------------	--------	----------------------	------------------------------------	------------------------------------	-----------------	----------------

PASS-THRU REFRIGERATORS

DH4-25PT	25	707	2°C to 8°C	23" x 30" x 59"	26" x 37.5" x 81.5"	CYCLE	450 LBS
DH4-50PT	50	1415	2°C to 8°C	48" x 30" x 59"	52" x 37.5" x 81.5"	CYCLE	600 LBS
DH4-80PT	80	2265	2°C to 8°C	75" x 30" x 59"	78" x 37.5" x 81.5"	CYCLE	800 LBS

STABILITY CHAMBERS

SCN4-24	23	651	4°C to 70°C	22" x 22.5" x 59"	27" x 26.5" x 81.5"	CYCLE	330 LBS
SCN4-24GD	23	651	20°C to 60°C	22" x 22.5" x 59"	27" x 36.5" x 81.5"	CYCLE	340 LBS
SCN4-52	49	1387	4°C to 70°C	49" x 25.5" x 59"	54" x 36" x 81.5"	CYCLE	510 LBS
SCN4-52GD	49	1387	20°C to 60°C	49" x 25.5" x 59"	54" x 36" x 81.5"	CYCLE	520 LBS

HUMIDITY STABILITY CHAMBERS

HCN4-24	23	651	4°C to 70°C	22" x 25.5" x 59"	26.75" x 35.75" x 80.75"	CYCLE	410 LBS
HCN4-25GD	23	651	20°C to 60°C	22" x 25.5" x 59"	26.75" x 35.75" x 80.75"	CYCLE	410 LBS
HCN4-52	49	1387	4°C to 70°C	49" x 25.5" x 59"	54" x 36.5" x 80.75"	CYCLE	610 LBS
HCN4-52GD	49	1387	20°C to 60°C	49" x 25.5" x 59"	54" x 36.5" x 80.75"	CYCLE	610 LBS