

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



LABORATORY RESEARCH & TESTING | PHARMACY STORAGE | MULTI-TEMPERATURE USE

Demanding Dual Storage

Our complete line of laboratory and pharmacy combination refrigerator/freezers provide the reliability and security essential for the protection of sensitive lab samples in the 0°C to -25°C range. Our scientific models offer stainless steel interiors and exteriors.

Seperated Cabinet Design

Refrigerator/Freezer combination units were designed specifically for use in a laboratory and pharmacy setting, offering many storage capabilities such as shelves and racks to store samples and other product with easy access to both freezer and refrigerator temperature ranges. Each compartment has it's own temperature control to regulate temperature.

Additional Security Options

So-Low Combination units also come with various options, such as dataloggers, chart recorders and audible visual alarm systems for temperature security.

Secure Protection

Units come standard with a door lock, protecting your valuable samples from tampering. Scientific models are UL Listed, and meet CDC recommendations for vaccine storage.

MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
-----------------	---------------	--------	----------------------	------------------------------------	------------------------------------	-----------------	----------------

GENERAL PURPOSE COMBINATION

DHK-10RFC (FRIDGE)	6.9	195	2°C to 8°C	23" x 20" x 26"	25.5" x 25.75" x 72.75"	AUTOMATIC	250 LBS
(FREEZER)	3.7	105	-10°C to -26°C	19" x 18" x 19"		MANUAL	
DHK-20RFC (FRIDGE)	11.8	335	2°C to 8°C	29.5" x 24" x 29"	36" x 29.25" x 77.75"	AUTOMATIC	450 LBS
(FREEZER)	7.2	205	-10°C to -26°C	25.5" x 22" x 22.25"		MANUAL	

SCIENTIFIC COMBINATION

DH-20RF (FRIDGE)	9.1	257	2°C to 8°C	28" x 28" x 20"	27.5" x 34" x 81.5"	AUTOMATIC	500 LBS
(FREEZER)	9.1	257	0°C to -23°C	28" x 28" x 20"		AUTOMATIC	
DH-48RF (FRIDGE)	21.01	595	2°C to 8°C	23" x 27" x 55"	52" x 33.75" x 81"	AUTOMATIC	575 LBS
(FREEZER)	21.01	595	0°C to -23°C	23" x 27" x 55"		AUTOMATIC	
DH-70RF (FRIDGE)	19.62	555	2°C to 8°C	22" x 28.5" x 55"	75" x 33.75" x 81.5"	AUTOMATIC	900 LBS
(FREEZER)	44.85	1270	0°C to -23°C	47" x 28.5" x 55"		AUTOMATIC	