

SO-LOW

ENVIRONMENTAL EQUIPMENT COMPANY

SETTING TEMPERATURE STANDARDS IN EVERY LAB

EST. 1959



ULTRA LOW FREEZER & LAB REFRIGERATION EQUIPMENT PRODUCT CATALOG

MEDICAL RESEARCH | INDUSTRIAL USE | COLD STORAGE | LAB TESTING | LIFE SCIENCES



LABORATORY RESEARCH & TESTING | CHEMICAL STORAGE | SECURE FMS PROTECTION

Spark Free Interior

Designed exclusively to control flammable material, these units do not contain any electrical components in the interior of the unit. The compressor electrical components have been sealed in a vapor-proof enclosure for additional safety.

NFPA & OSHA

Our Flammable Refrigerators and Freezers are designed to meet National Fire Protection Agency NFPA and OSHA guidelines for the refrigerated storage of flammable (volatile) material storage in general laboratory work areas.
Meets NFPA and OSHA guidelines for hazardous locations, Class 1, Division I, Division II, Group C & D.

Electrical Listing

Units are UL or ETL listed. Varies by model.

MODEL NUMBER	CUBIC FEET	LITERS	TEMPERATURE RANGE	INTERIOR DIMENSIONS W x F-B x H	EXTERIOR DIMENSIONS W x F-B x H	DEFROST TYPE	UNIT WEIGHT
-----------------	---------------	--------	----------------------	------------------------------------	------------------------------------	-----------------	----------------

REFRIGERATOR

MV4-6UCFMSR	5	141	1°C to 10°C	20.25" x 17" x 27"	23.75" x 24.25" x 32.25"	MANUAL	75 LBS
DHH4-14SDFMS	14	396	1°C to 10°C	19" x 16.75" x 63.75"	23.5" x 24.25" x 73"	MANUAL	150 LBS
DHH4-17SDFMS	17	481	1°C to 10°C	26.25" x 20.75" x 58.5"	30.75" x 28.33" x 67.75"	MANUAL	190 LBS
DHH4-20SDFMS	20	566	1°C to 10°C	24.75" x 20.75" x 65.25"	30.75" x 29" x 74.5"	MANUAL	210 LBS
DHH4-30SDFMS	30	849	1°C to 10°C	32" x 22.5" x 65"	35.75" x 35.25" x 72.25"	MANUAL	330 LBS

FREEZER

DHH20-14SDFMS	14	396	-15°C to -25°C	19" x 17.5" x 63.75"	23.75" x 24.75" x 73"	MANUAL	145 LBS
DHH20-17SDFMS	17	481	-15°C to -25°C	26.25" x 22.5" x 58.25"	30.75" x 28.5" x 67.75"	MANUAL	170 LBS
DHH20-20SDFMS	20	566	-15°C to -25°C	25.75" x 21.75" x 65"	30.75" x 29" x 74.5"	MANUAL	180 LBS
DHH20-30SDFMS	30	849	-18°C to -25°C	32" x 22.5" x 65"	35.75" x 35.5" x 72.5"	MANUAL	330 LBS

REFRIGERATOR / FREEZER COMBINATION

DHH-9RFFMS	(REF)	4.6	130	1°C to 10°C	20.25" x 17" x 27"	23.75" x 24.25" x 65.5"	MANUAL	150 LBS
	(FRZ)	4.4	125	-15°C to -25°C	19" x 16" x 26"			
DHH-30RFFMS	(REF)	25	707	1°C to 10°C	29" x 22.25" x 12"	35.75" x 35.5" x 72.5"	MANUAL	350 LBS
	(FRZ)	5	141	-15°C to -25°C	32" x 22.5" x 50.25"			