

Access® High-Low PLUS
Power Exam Table
with 17" Low Seat Height

Setting the Bar for Accessibility



Meet the ADA-compliant Access PLUS with 17" low transfer height – designed for today's standards and tomorrow's recommendations. Stay compliant and future-ready.



Discover how Access PLUS
meets the standard.



Create an elevated patient-provider experience.



**2024 Access Board ADA
Recommended 17" Seat Height**




**ADA-Compliant ComfortGrip
Assist Handles**



**ADA-Compliant Two-Side
Transfer Access**



**ADA-Compliant Seat Transfer
Surface - 27.5" W x 17" D**



Up to 700 lb. Weight Capacity



**Industry-Best 6 Cubic Feet
of Storage**



**Women's Health Package: Pelvic
Tilt, Stirrups & Drawer Warmer**



**Whisper-Quiet, Effortless
EasyGlide Mobility System**



**Generous Front Storage Drawer
with Heater Option**



The Brewer Difference.

Industry's most
safety features.

Industry's only
5-year warranty.

USA-made in
the Midwest.

Access® High-Low PLUS Power Exam Table with 17" Low Seat Height

Specifications

Advanced Safety

- Exclusive SafeGlide™ leg extension increases safety and satisfaction
- ComfortGrip™ Patient Assist Handles facilitate easy transfers and adjustments, enhancing patient independence and safety during examinations
- Return to chair feature on power back models
- Anti-tamper on/off switch limits accidental motion
- Industry-best three safety zones offer protection against pinching and impact-related injuries.
- 17" (43.2 cm) low height that raises to 37" (94 cm), both ADA compliant
- Easy on/off for patients
- Powered backrest adjusts for easy patient positioning*

** Pneumatic back models are also available*

Patient Comfort

- Easy lateral wheelchair transfers
- Multi-density support upholstery
- Toe/leg extension & stainless steel debris tray

Clinician Convenience

- Signature Storage Hub with patented pass-through side drawers that are accessible from right or left side of the table
- Flip-top storage compartment completes the storage hub offering additional space to keep paper rolls, supplies and equipment.
- **Front drawer:**
19.75" x 19.25" x 4" (50.2 cm x 48.9 cm x 10.2 cm)
- **Side drawer (each):**
22.5" x 20.25" x 8" (57.2 cm x 51.4 cm x 20.3 cm)
- **Paper roll storage capacity:**
23" x 23.25" x 3.25" (58.4 cm x 59.1 cm x 8.3 cm)
- **Slide-through paper roll holder:** accepts 18" (45.7 cm) and 21" (53.3 cm) rolls

Table Specifications

- 450 lb (204 kg) weight capacity
(6000PLUS-17, 6001PLUS-17, 6500PLUS-17, 6501PLUS-17)
- 700 lb (317 kg) weight capacity
(6800PLUS-17, 6801PLUS-17)
- ComfortGrip™ patient assist handles come standard
- Ergonomic adjustable stirrups come standard
- Foot control comes standard
- Industry best 5-year warranty
- **Height:** 17"-37" (43.2 cm - 94 cm) — Both heights meet the 2024 Access Board recommendation.
- **Length:** 60", extends to 78" (152.4 cm - 198.1 cm)
- **Upholstery Width:** 27.5" (69.9 cm)
- **SafeGlide™ Leg extension length:** 18" (45.7 cm)
- 115VAC, 60Hz outlet
- Safe, low voltage actuators
- Up to 6° pelvic tilt
- 60-Watt front drawer heater
- Easy to clean, seamless upholstery

Certifications

All Brewer products conform to US & Canadian electrical & medical industry codes & regulations and meet or exceed industry requirements.

- UL 60601-1
- CAN-CSA C22.2 No. 601.1
- ISO 13485:2016



17" Seat Height Option:

- **6000PLUS-17**
Pneumatic back with stirrups, assist handles & 17" low-transfer height
- **6001PLUS-17**
Pneumatic back with stirrups, assist handles, pelvic tilt, drawer warmer & 17" low-transfer height
- **6500PLUS-17**
Power back with stirrups, assist handles & 17" low-transfer height
- **6501PLUS-17**
Power Back w/outlet, stirrups, assist handles, pelvic tilt, drawer warmer & 17" low-transfer height
- **6800PLUS-17**
Power back with stirrups, assist handles & 17" low-transfer height
- **6801PLUS-17**
Power Back w/outlet, stirrups, assist handles, pelvic tilt, drawer warmer & 17" low-transfer height

Accessories:

- **3500687PLUS:**
EasyGlide™ Mobility System
- **99502:**
One Universal Socket Clamp
- **99504:**
Deluxe Articulating Knee Crutches*
- **100253-TG:**
Urology Drain Pan System
- **100407PLUS:**
Front Drawer Dividers
- **100540:**
Two Support Adapters for Knee Crutches
(factory installed)
- **100550:**
Two Support Adapters for Knee Crutches
(field installed)
- **102031PLUS:**
Brewer Welch Allyn® IV Bracket
(factory installed)
- **102032PLUS:**
Brewer Welch Allyn® IV Bracket
(field installed)
- **101800PLUS:**
6001PLUS-17 Outlet

**Requires two 99502 and one of either 100540 or 100550*



Create an elevated patient-provider experience.

Models	6000PLUS-17	6001PLUS-17	6500PLUS-17	6501PLUS-17	6800PLUS-17	6801PLUS-17
Weight Capacity	450	450	450	450	700	700
Access Board 17" Seat Height Recommendation	✓	✓	✓	✓	✓	✓
Power Base	✓	✓	✓	✓	✓	✓
Power Back			✓	✓	✓	✓
Manual Back	✓	✓				
Patient Assist Handles	✓	✓	✓	✓	✓	✓
Stirrups	✓	✓	✓	✓	✓	✓
Pelvic Tilt & Drawer Warmer		✓		✓		✓
Outlet		Optional		✓		✓
EasyGlide Mobility	Optional	Optional	Optional	Optional	Optional	Optional
Mobility	Optional	Optional	Optional	Optional	Optional	Optional

Available Colors.

 Chateau Grey AF906	 Flannel AF945	 Horizon AF928	 Mist AF938	 Black Satin PR40
 Cool Grey MV821	 Eclipse AF919	 Blue Fog PR46	 Cocoa PR51	 Cabernet PR36
 Feather PR37	 Deep Sapphire PR55	 Gunmetal PR92	 Clamshell PR57	 Ivy PR53



The Brewer Difference.

Industry's most
safety features.

Industry's only
5-year warranty.

USA-made in
the Midwest.