

HUNTLEIGH

A MEMBER OF THE ARJO FAMILY



Dopplex

Vascular Assessment
Accessories & Consumables
Catalogue

2023



CONTENTS

<u>DOPPLER ACCESSORIES</u>	<u>A</u>
High Sensitivity Probes	3
Intraoperative Probes / Packs	4
Generic Accessories	5
Generic Consumables	9
 <u>ABILITY ACCESSORIES</u>	 <u>B</u>
Cuff Kits	10
Consumables	12
 <u>LEGACY ACCESSORIES</u>	 <u>C</u>
Dopplers	13

VP4XS

High Sensitivity 4MHz probe
(Detection of Deep Lying vessels)

**VP5XS**

High Sensitivity 5MHz probe
(Oedematous Limbs & Deep Lying Vessels)

**VP8XS**

High Sensitivity 8MHz probe
(Detection of Peripheral Vessels & Calcified Arteries)

**VP10XS**

High Sensitivity 10MHz probe
(Detection of Smaller Vessels)

**EZ8XS**

High Sensitivity Wide-beam 8MHz probe
(Detection of Peripheral Vessels & Calcified Arteries)

**PPGA1**

Arterial PPG Adaptor and Probe
(Detection of Toe Pressures - Used in ATP Kits)



Compatible



Non-Compatible

DIPP10

Pre-sterilised Single-Use Intraoperative 8MHz Doppler probes (Box of 10)

DMX	SD2	D900
☒	☒	☒



DISP10XS

Single-Use Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Pole Clamp & DIPP10

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



DISP10-DMXR

Single-Use Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Pole Clamp & Box of 10 Single-use Intraoperative 8MHz Doppler Probes and DMXR Doppler

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



PA8XS

Intraoperative Probe Adapter - DMX range

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ISP-XS

Intraoperative Probe Starter Pack: Includes PA8XS Probe Adaptor, Clamp

☒	☐	☐
-------------------------------------	--------------------------	--------------------------



ACC47

IV Pole Clamp for use with Intraoperative Probe

☒	☒	☒
-------------------------------------	-------------------------------------	-------------------------------------



ACC21

Stereo Headset

**ACC34**

Soft Carry Pouch

**ACC232**Dopplex Stand - shipped flat packed,
Includes Adaptor Plate**ACC159**Adult Arm/Ankle Cuff with Colder connector
(23-33cm)**ACC160**Large Arm/Ankle Cuff with Colder connector
(31-40cm)**ACC161**Thigh Cuff with Colder Connector
(38-50cm)

Compatible



Non-Compatible

ACC162

Large Digit/Toe Cuff with Colder Connector
(2.5cm) (pack of 5)

ACC163

Standard Digit/Toe Cuff with Colder Connector
(1.9cm) (pack of 5)

ACC185

Digit /Toe Cuff with colder connector
- 3 Large & 3 Standard

ACC164

Tourniquet Cuff with Colder Connector

ACC168

Aneroid Sphygmomanometer with trigger
release valve and Colder Connector

ACC187

Tourniquet cuff & Aneroid Sphygmomanometer
Kit with trigger release valve & Colder Connector

ACC188

Blood Pressure cuff & Aneroid
Sphygmomanometer Kit with trigger release
valve & Colder Connector

DMX
SD2
D900



ACC226

Medical Grade Rechargeable Kit

DMX
SD2
D900**ACC227**

Micro SD card

**ACC228**

Doppler Stand Adaptor plate

**ACC-VAS-028**

T-Tube Toe Cuff Inflator System

**ACC-VAS-029**

Large Carry bag - use with ATP/ABI Kits



Compatible



Non-Compatible

ACC-VAS-030

USB cable use with Vascular Reporter (3m)

DMX	SD2	D900
☒	☐	☐



ACC-OBS-080

Splash-proof Doppler pouch (IPX2)

DMX	SD2	D900
☒	☐	☐



ACC3

Ultrasound Gel
(Box 12 x 250ml Bottles)

DMX	SD2	D900
☒	☒	☒



ACC24

Ultrasound Gel
(Box 12 x 60ml Tubes)

DMX	SD2	D900
☒	☒	☒



ACC-1300-0145

Ultrasound Gel (box 48 x 20gm Sachet)

DMX	SD2	D900
☒	☒	☒



APPG Adhesive Straps
(Pack of 240)

Neuropen & Tip
- Includes 10g Monofilament

	<i>DMX</i>	<i>SD2</i>	<i>D900</i>



DOPPLEX ABILITY



ACC-VAS-007

Large Adult Right Arm Cuff - Red tube
(34-46cm)

RIGHT	LEFT	KIT
☒	☐	☐



ACC-VAS-008

Large Adult Right Ankle Cuff - Black tube
(24-35cm)

RIGHT	LEFT	KIT
☒	☐	☐



ACC-VAS-009

Large Adult Left Arm Cuff - Yellow tube
(34-46cm)

RIGHT	LEFT	KIT
☐	☒	☐



ACC-VAS-010

Large Adult Left Ankle Cuff - Green tube
(24-35cm)

RIGHT	LEFT	KIT
☐	☒	☐



B

☒ Compatible

☐ Non-Compatible

ACC-VAS-011

Large Adult Cuff Set
(34-46cm)

RIGHT	LEFT	KIT
☐	☐	☒

**ACC-VAS-023**

Adult Right Arm Cuff - Red tube
(22-36cm)

RIGHT	LEFT	KIT
☒	☐	☐

**ACC-VAS-024**

Adult Right Ankle Cuff - Black tube
(18-28cm)

RIGHT	LEFT	KIT
☒	☐	☐

**ACC-VAS-025**

Adult Left Arm Cuff - Yellow tube
(22-36cm)

RIGHT	LEFT	KIT
☐	☒	☐

**ACC-VAS-026**

Adult Left Ankle Cuff - Green tube
(18-28cm)

RIGHT	LEFT	KIT
☐	☒	☐

**ACC-VAS-027**

Adult Arm Cuff Set
(22-36cm)

RIGHT	LEFT	KIT
☐	☐	☒



Compatible



Non-Compatible

ACC-VAS-012

Fixing Plate (required for ACC-VAS-013 or ACC-VAS-154)



ACC-VAS-013

Pole Stand / Trolley
(Requires ACC-VAS-012)



ACC-VAS-015

Carry bag



ACC-VAS-154

Wall Mount (requires ACC-VAS-012)



ACC-VAS-016

Infection Control Sleeves
(Box of 100)



ACC-VAS-017

Standard Thermal Paper - Pack of 5 (approx. 55 prints per roll)



ACC-VAS-019

Self-adhesive Thermal Paper - Pack of 5
(approx. 55 prints per roll)



LEGACY



ACC171

Rechargeable Battery Kit (UK)

ACC172

Rechargeable Battery Kit (USA)

ACC173

Rechargeable Battery Kit (EU)

-Supplied with Charger & 2 Batteries

D900
SD2
MD200



PA8-HG

Intraoperative Probe Adapter



DISP10

Single-Use Intraoperative Probe Starter Pack:
Includes PA8-HG Probe Adaptor, Pole Clamp
& Box of 10 Single-use Intraoperative 8MHz
Doppler Probes



ISP

Intraoperative Probe Starter Pack:
Includes PA8-HG Probe Adaptor, Clamp



Compatible



Non-Compatible

D900
SD2
MD200



ACC232

Dopplex Stand - Shipped flat packed

ACC180

Dopplex Stand - Shipped flat packed

ACC190

RS232-USB Adaptor for use with Dopplex Reporter Software Package (DR4), running on a PC without a RS232 port - requires ACC35 Cable