

When a Rollator becomes a Hybrid...



Introducing the Lumex® HybridLX Rollator Transport Chair

The HybridLX is the latest addition to the Lumex line of Mobility products, a line that is well known and respected for its innovation and quality.

The HybridLX is an innovative rollator that combines the function of a wheeled walker and a transport chair into one lightweight design. Users, when they need a rest, can easily transform this rollator into a transport chair with easy mount swingaway footrests. Keeping the patient's desire for independence and comfort in mind, the HybridLX is designed with secure handgrips, an easy braking system and a wide seat and backrest for comfort.

The HybridLX is ideal for improving mobility and maintaining an active and independent lifestyle.

FEATURES:

- **Hybrid concept** that easily transforms the rollator into a transport chair with removable footrests (adjustable in length and with heel loops) when user needs a short rest
- **Lightweight aluminum frame** that is durable weighs only 22 lbs (this includes footrests) making it easier for users to move about freely
- **Ultra-Compact design** is easy to assemble without the need for tools and folds effortlessly for convenient transport and storage in the optional carry bag
- **Easy-to-use handbrakes** allow for short and longer period braking when used as rollator and transport chair (user must be accompanied when using as a transport chair)
- **Backrest and 18" wide seat** are flexible and made of soft polyester with pull up seat strap for easy folding provide added comfort for user
- **Large 8" wheels** for easy rolling indoor and outdoor on rough surfaces
- **Storage pouch** included for personal items
- **Attractive, contemporary colors** available: Titanium (Item LX1000T) and Majestic Blue (Item LX1000B)



graham field™

The Graham-Field® Family of Brands

LUMEX Everest & Jennings **10101** Hausted® **Grafo** **Basic American** **LUMISCOPE** **JOHN BUNN** **SIMMONS**

SIMMONS Clinical Care products manufactured by GF Health Products, Inc. under license from the trademark owner, Deemedwell, LLC, a subsidiary of Simmons Bedding Company, Atlanta, Georgia. SIMMONS is a trademark of Deemedwell, LLC, a subsidiary of Simmons Bedding Company. © 2010 Deemedwell, LLC. All rights reserved.