

Stainless Steel Deluxe IV Stand

- Constructed of chrome-plated stainless steel
- 5 leg IV stand with 2" nylon swivel casters provides easy movement in various patient care settings
- Telescoping pole adjusts from 51" – 86"
- 4 circular hooks for solution suspension
- 2 year limited warranty
- Designed to meet the requirements of HCPCS code: E0776

7016A

1 ea

Specifications

Height Adjustment:51" - 86"

Product Weight:7.5 lb

Warranty:2 years



GRAHAM-FIELD