

ROLLATORS

HYBRID LX ROLLATOR TRANSPORT CHAIR

The HybridLX is an innovative rollator that combines the function of a wheeled walker and a transport chair into one lightweight design. Keeping the patient's desire for independence and comfort in mind, the HybridLX is designed with secure handgrips, a comfortable and wide seat and backrest, and the effortless design to transition from a rollator to a transport chair.

Features

- **Lightweight frame** – durable, aluminum frame weighs only 22 lbs (this includes footrests) making it easier for users to move about freely
- **Hybrid concept** – designed to be used primarily as a rollator while easily transforming into a transport chair with removable footrests (adjustable in length) when user needs a short rest
- **Ultra-Compact design** – easy to assemble without the need for tools and folds effortlessly for convenient transport and storage
- **Easy-to-use handbrakes** – allow for short and longer period braking when used as rollator and transport chair (user must be accompanied when using as a transport chair)
- Flexible, soft polyester backrest and 18" wide mesh seat with pull up strap for easy folding provides comfort for user
- Large 8" wheels for easy rolling indoor and outdoor on rough surfaces
- Complete with storage pouch for personal items
- Available in attractive, contemporary colors: Titanium and Majestic Blue
- Designed to meet the requirements of HCPCS Codes: E0143 + E0156
- Three (3) year limited warranty



		overall Width		overall Depth with footrests		Depth front to rear		Width inside handgrips		minimum Handle Height		maximum Handle Height		Wheel diameter	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
LX1000T	Titanium	26	66.04	40	101.60	28.5	72.39	18.5	46.99	30	76.20	37	93.98	8	20.32
LX1000B	Majestic Blue	26	66.04	40	101.60	28.5	72.39	18.5	46.99	30	76.20	37	93.98	8	20.32

		Seat Width		Seat Depth		Seat Height from floor		user Height Range		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	feet	m	lb	Kg	lb	Kg	
LX1000T	Titanium	18	45.72	11	27.94	21.25	53.98	5' 4" - 6' 2"	1.63 - 1.88	22	9.98	300	136.08	1 ea
LX1000B	Majestic Blue	18	45.72	11	27.94	21.25	53.98	5' 4" - 6' 2"	1.63 - 1.88	22	9.98	300	136.08	1 ea



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

© 2008, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.