

## Deluxe Tilt Overbed Table

- Chrome plated and welded tubular steel column and base make the table durable and stable.
- 1 ½" hooded, ball-bearing casters provide easy maneuverability across all floor surfaces.
- Table top adjusts 13" in infinite positions and is constructed of wood-grained oak Melamine laminate on particle-board core.
- Raised T-molded plastic edge.
- Chrome-plated "H" base and column for aesthetics and durability.
- 2 locking casters
- Designed to meet the requirements of HCPCS code: E0274.
- Supports up to 30 lb of dispersed weight.
- Warranty is Three years

**GF8905-1A**      1 ea



| Item      | Table Height | Table Width | Table Depth | Base Width | Base Depth | Weight | Dispersed Weight Capacity |
|-----------|--------------|-------------|-------------|------------|------------|--------|---------------------------|
| GF8905-1A | 30" – 43.8"  | 30"         | 15.2"       | 26.2       | 16.5       | 21 lb  | 30 lb                     |