

CLEANING & CARE GUIDELINES

INSTRUCTIONS OF USE

COPYRIGHT ©2025. ALL RIGHTS RESERVED



USE

Phillips Safety Products must be used only for their designated purposes.

WARRANTY

Phillips Safety Radiation Apparel products are covered by a two-year warranty against manufacturing defects.

Attention: This does NOT include normal wear and overloading of use.

CLEANING

Use a mild detergent or cleaner to remove any dirt or debris.

Do NOT use abrasives.

DISINFECTING

It can be disinfected with most standard hospital-grade disinfectants containing QATs, IPA, or Hydrogen Peroxide.