

REF 6001

Buffing Pad

Intended Use

The Buffing Pad, designed for use on rigid endoscopes, features a buffing surface that helps reduce debris from the lens.

Warnings

- Use only with rigid endoscopes.
- DO NOT use on very hot surfaces.
- DO NOT use with chemicals or solvents.
- DO NOT make any modifications to pads.
- DO NOT bend scope while buffing

In The Event Of An Incident

Any serious incident occurring in connection with this device must be reported to the manufacturer and the relevant authority of the member state in which the user and/or patient is based.



