

PRODUCT DATASHEET

Model: VFS

Description: DR Weight Bearing Foot Stool

Revision: November 2020

Product Attributes

- **For standing foot exams**
- **Handrails for patient support**
- **Single step platform**
- **Lightweight**
- **To be used with a U-Arm**



DIMENSIONS

Base Dimensions: 35x14x24 in. (88.9 x 35.56 x 60.96cm)

Handrail height: 35 in. (88.9 cm)

Weight: 33.6 lbs. (15.24 kg)

Weight Tolerance: 500 lbs. (226.8 kg)

U-Arm Cavity: 10.25 x 28.5 in. (26 x 72.4 cm)

Imaging Area: 17.25 x 17.25 in. (43.8 x 43.8 cm)

Boxed Dimensions and Weight

Box 1: 60 x 26 x 15" (152.4 x 66 x 38.1 cm), Boxed Weight: 75 lbs. (34 kg)

FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

Country	Class
USA	FDA Class 1, 510(k) exempt
CANADA	Class 1
EUROPE	Class 1, Annex VII, CE
TGA	Class 1