

PRODUCT DATASHEET

Model: POSPLC

Description: Positioning Partner Lite

Revision: January 2026

Product Attributes

- Counterbalanced Vertical travel
- Small footprint
- Easy to drape
- No electric
- Locking casters
- Simple one-handed height adjustment

Dimensions

Product Dimension: 28.3 x 18.8 x 58.3 in. (71.8 x 47.8 x 148.1 cm)

Total Vertical Stroke: 39 in. (99.06 cm)

Weight: 45 lbs. (20.4kg)

Boxed Dimensions: 24 x 20 x 57 in. (60.96 x 50.8 x 144.78 cm)

Boxed Weight: 75 lbs. (34kg)

Flammability Rating

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

Biocompatibility

This model is a biocompatible device per ISO 10993.

Statement Of Conformance

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

Registrations and Compliance

The model listed above is made in the USA and registered to:

Reina Imaging



*Protect-A-Grid Encasement not included