

PRODUCT DATASHEET

Model: **SPH**

Description: **Stitching Platform with Handles**

Revision: June 2015

PRODUCT ATTRIBUTES:

- Patient immobilizer
- Handrails for stability
- Slip resistant platform
- Casters

DIMENSIONS AND PERFORMANCE


Dimensions: 21 x 26.125 x 75.25 in. (53.34 x 66.36 x 191.14 cm)

Base Height: 9 in. (22.86 cm)

Inside Width: 24 in. (60.96 cm)

Weight: 61 lbs. (27.67 kg)

Weight Tolerance: 750 lbs. (340 kg)



FLAMMABILITY RATING

The product is non-flammable under the normal operating conditions of the DR Panel.

UL Rating: UL94V-2.

BIOCOMPATIBILITY

This model is a biocompatible device per ISO 10993.

STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

DEVICE FAMILY CLASSIFICATION

Country	Class
USA	FDA Class 1, 510(k) exempt
Canada	Class 1
Europe	Class 1, Annex VII, CE
TGA	Class 1