

WHEELCHAIRS



Specialty Folding Wheelchairs

Advantage® Recliner

3K01 & 3K02 Series

The Advantage Recliner, with a higher back and removable, full-support headrest for additional neck support, adjusts easily to fit every need.

- Available in 16, 18 and 20 inch (40.64, 45.72 and 50.80 cm) seat widths, 17 inch (43.18 cm) seat depth and in 22 inch (55.88 cm) seat width, 19 inch (55.88 cm) seat depth
- Durable hammertone frame features manual full reclining back range from 90-170 degrees
- Comfortable desk and full length arms with padded armrests
- Metal spreader bar on back which works as anti-fold device
- Embossed black vinyl upholstery with standard chart pocket on back
- Anti-tippers standard
- Tool-free adjustable elevating front legrests
- 450 lb (204.12 kg) weight capacity on 22 inch wide models - EVENLY DISTRIBUTED
- 300 lb (136.08 kg) weight capacity on 16, 18 & 20 inch wide models - EVENLY DISTRIBUTED
- Five (5) year limited warranty



Rehab Shower Commode

12022 & 12023 Series

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. Designed with a corrosion-resistant stainless steel frame, the chair can be positioned over a standard toilet or used as a portable, self-contained commode.

- Available with durable 24 inch (60.96 cm) rear wheels with handrims and solid rubber tires
- Removable full length arms with composite armrests
- Available in Low-back 14 ¼ inch (36.83 cm) on 16 inch (40.64 cm) and Hi-back 17 ¼ inch (43.82 cm) on 19 ½ inch (48.26 cm)
- Available in 17 inch (43.18 cm) seat with 18 inch (45.72 cm) seat depth on Low-back 14 ¼ and 17 ¼ inch and the wider version 20 inch (50.80 cm) seat width on 16 inch Low-back and Hi-back 19 ½ inch models
- Sealed vinyl cushioned oval seat and frame side contours allow easy access and comfort in toileting activities and bowel management procedures
- Removable, hygienic plastic pail and lid are included
- Non-corrosive composite footplates
- 250 lb (113.40 kg) weight capacity -EVENLY DISTRIBUTED on 14 ¼ and 17 ¼ inch back models
- 300 lb (136.08 kg) weight capacity -EVENLY DISTRIBUTED on 16 and 19 ½ inch back models
- Two (2) year limited warranty

