




Advantage[®]

Recliner

PARTS CATALOG



CONTENTS

FRAME AND CROSSBRACE 3

ARM 4

CASTER..... 4

UPHOLSTERY 5


FRONT RIGGING, DETACHABLE 6

REAR WHEEL 7

ACCESSORIES 7

Returns: Parts may not be returned without prior authorization from GF Health Products, Inc. Custom parts cannot be returned for credit. Please refer to GF Health Products, Inc. [Terms and Conditions](#) for further information.

FRAME AND CROSSBRACE



1 Sideframe

16" x 17", 18" x 17", 20" x 17"

Right..... 90763701

Left 90763702

22" x 19"

Right..... 90763707

Left 90763708

2 Crossbrace, pair; includes center bolt and nut

16" wide x 17" deep 907637A1

18" wide x 17" deep 907637A2

20" wide x 17" deep 907637A3

22" wide x 19" deep 907637A4

3 Bottom inner rail kit; includes bottom inner rail, 2 each; mounting screw and nut, 4 each; and sleeve, 2 each..... 90763703

4 Plug, seat rail, set of 4 90763704

5 Grip, black, 10 each 90763005

6 Anti-tipper, pair

16" x 17", 18" x 17", 20" x 17" 90763730

22" W x 19" 90763067

7 Seat slide, front and rear

16" - 18" wide, pair 90763705

20" wide, set 90763706

22" wide, set 90763709

8 Backpost, cable and trigger, pair; includes cable mounting hardware 90763710

9 Mech-lock, pair; includes mounting hardware..... 90763715

10 Spreader bar knob, pair..... 90763725

11 Spreader bar

16" wide 907637E1


18" wide 907637E2

20" wide 907637E3

22" wide 907637E4


12 Headrest post, pair..... 90763720

ARM




- 1 Arm, detachable fixed-height, pair:
Includes right and left arm with side
panel and spring lock (does not
include arm pad)
Desk length..... 90763020
Full length..... 90763021
- 2 Arm pad, pair: Includes
two black arm pads
and mounting screws
Desk length..... 90763022
Full length..... 90763023
- 3 Spring lock (10 each) 90763024

CASTER



- Complete caster (fork and wheel),
8" x 1" solid rubber tire
- 16" x 17", 18" x 17", 20" x 17" 90763560
22" W x 19" 90763570
- 1 Caster wheel, axle and nut
16" x 17", 18" x 17", 20" x 17" 90763561
22" W x 19" 90763571
 - 2 Caster fork, stem nut and dust cap
16" x 17", 18" x 17", 20" x 17" 90763562
22" W x 19" 90763572
 - 3 Caster stem bearing (10 each)
16" x 17", 18" x 17", 20" x 17" 90763563
22" W x 19" 90763573
 - 4 Caster dust cap (10 each)
A 16" x 17", 18" x 17", 20" x 17" 90773544
B 22" W x 19" 90763044
 - 5 Caster axle and nut (10 each)
16" x 17", 18" x 17", 20" x 17" 90763564
22" W x 19" 90763574



1 Back upholstery, black vinyl, with built-in chart holder: Includes hardware

16" wide	907637B1
18" wide	907637B2
20" wide	907637B3
22" wide	907637B4

2 Seat upholstery, black vinyl: Includes hardware

16" wide x 17" deep	907637C1
18" wide x 17" deep	907637C2
20" wide x 17" deep	907637C3
22" wide x 19" deep	907637C4


3 Headrest upholstery, black vinyl: Includes hardware

16" wide	907637D1
18" wide	907637D2
20" wide	907637D3
22" wide	907637D4

4 Headrest cushion, black vinyl..... 90763735


5 Carry strap, black vinyl, pair: Includes hardware..... 90763740

FRONT RIGGING, DETACHABLE




Footrest, pair, all widths
Includes right and left hanger with release
and composite footplate with extension
and heel loop 90763030

Composite footplate, pair
(Fits all front rigging) 90763033




Low-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
composite footplate with extension, and
legrest panel (does not include
heel loop) 90763031

Elevating legrest panel kit
(Fits all elevating legrests)
Includes 2 black padded panels
and hardware 90763032



High-pivot elevating legrest, pair, all widths
Includes right and left hanger with release,
aluminum footplate with extension, and
legrest panel (does not include
heel loop) 90763435

REAR WHEEL



1 Rear wheel with handrim, 24" x 1", solid rubber tire: Includes hub cap and bearing; does not include axle or nut	90763550
2 Hub cap (10 each).....	90763051
3 Wheel bearing (10 each)	90763552
4 Rear axle and nut (10 each).....	90763554

ACCESSORIES

Advantage Recliner accessories (not shown)

1 Heel loop, pair.....	90763064
2 Wheel lock, pair: Includes rubber tips	
Push to lock.....	90763060
Pull to lock.....	90763061
3 Rubber tip, wheel lock (10 each).....	90763063
4 Wheel lock extension, pair	90763062
5 Positioning belt	
Airline-style metal buckle, adjusts 27"-47"	92133040
Auto-style push-button, adjusts 27"-37"	300-505
6 Anti-Fold bar	92523850