

WHEELCHAIRS



Lightweight Folding Wheelchairs

Traveler® L3 Plus and Traveler® L3 Plus XVI

3F012, 3F013 & 3F014 Series

The Traveler L3 Plus is an improved version of the previous Traveler L3 series. The wheelchair comes standard with adjustable anti-tippers and a fully detachable 2 inch (5.08 cm) seat extension with hardware.

- Available in 16, 18 and 20 inch (40.64, 45.72 and 50.80 cm) seat widths including 2 inch (5.08 cm) removable seat extension kit converting the chair from 16 inch depth to 18 inch depth if necessary
- Two position axle and caster fork provides 18 to 20 inch seat height in regular models and seat height 16.30 to 17.90 inch (41.40 to 45.47 cm) in the XVI models
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Non-marking polyurethane tires with 24 inch (60.96 cm) rear wheels with 8 inch (20.32 cm) casters in the regular models and 20 inch rear wheels with 6 inch (15.24 cm) casters in the XVI models
- Anti-tippers come standard
- Flip-back padded arms available in fixed desk length, adjustable height desk and detachable flip-back full length
- Hammertone finish and carbon steel frame provide strength and durability
- Tool-free adjustable front rigging with heel loops standard on swingaway footrests
- 300 lb (136.08 kg) weight capacity - EVENLY DISTRIBUTED
- Five (5) year limited warranty



24 inch (60.96 cm) wheel model shown



24 inch (60.96 cm) wheel model shown



20 inch (50.80 cm) wheel model shown

Traveler® L4

3F02 & 3F03 Series

The Traveler L4 is a high-strength lightweight chair that weighs less than 34 lb (15.42 kg) [without front rigging] with an adjustable back.

- Available in 16, 18 and 20 inch (40.64, 45.72 and 50.80 cm) seat widths and 16 inch seat depth. Plus 18 and 20 inch seat widths, 18 inch seat depth
- Quick release axle and height adjustable desk arms are available for some models
- Black nylon upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Anti-tippers come standard
- Comfortable flip-back desk and full length arms with padded armrests
- 250 lb (113.40 kg) weight capacity - EVENLY DISTRIBUTED
- 300 lb (136.08 kg) weight capacity - EVENLY DISTRIBUTED on models with 18 inch seat depth
- Five (5) year limited warranty



K0004*