

WHEELCHAIRS





Everest & Jennings®

by graham-field

Value:

Graham-Field is known for a broad product offering at affordable pricing, without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell to and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings® line of wheelchairs.

Quality:

The E&J™ name has been synonymous with quality manual folding wheelchairs since 1933. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, Graham-Field is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Georgia; St. Louis, Missouri; and Los Angeles, California. Graham-Field also carries a full line of E&J parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.



Bariatric Folding Wheelchairs

Traveler® HD

3G01 Series

The Traveler HD is a heavy duty wheelchair featuring a reinforced frame, double crossbrace, and a long lasting, embossed back and seat vinyl upholstery with a double inner liner for durability and support.

- Available in 20, 22 and 24 inch (50.80, 55.88 and 60.96 cm) seat widths, 18 inch (45.72 cm) seat depth (1 inch / 2.54 cm longer than standard)
- Embossed black vinyl upholstery with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Reinforced frame at all critical load-bearing points
- Two-position axle and caster fork provides 2 inch (5.08 cm) seat height adjustment
- Tool-free adjustable front rigging with standard heel loops on footrests and durable aluminum black footplates
- Comfortable full and desk length arms with padded armrests
- 500 lb weight capacity - EVENLY DISTRIBUTED
- Five (5) year limited warranty



Traveler® HTC

3J01 Series

The Traveler HTC (Hospital Transport Chair) is an all-in-one solution wheelchair for the clinical transport need, as it incorporates a special attendant push handle and rear mounting attendant wheel locks for better caregiver assistance.

- Fixed desk arms and front rigging lock nuts ensure inventory control
- Anti-tippers and one O2/IV telescopic holder come standard
- 500 lb (226.80 kg) weight capacity -EVENLY DISTRIBUTED
- Five (5) year limited warranty



Paramount XD

5PX10 Series

The Paramount XD, built for the bariatric wheelchair user, is constructed of 1 inch (2.54 cm) tubing to meet high weight capacity demands.

- Available in 26 and 30 inch (66.04 and 76.20 cm) seat widths, 19 ½ inch (49.53 cm) seat depth
- Comfortable desk arms with padded armrests
- Front rigging with extra wide aluminum footplate
- 8x2 inch (20.32x5.08 cm) non-marking solid casters and 24 inch (60.96 cm) solid rear wheels with chromed handrims
- Black reinforced nylon upholstery with standard chart pocket on back
- 650 lb (294.84 kg) weight capacity - EVENLY DISTRIBUTED
- Five (5) year limited warranty

