

## BARIATRIC/HEAVY DUTY WHEELCHAIRS

# TRAVELER® HD WHEELCHAIR

The Traveler® HD is manufactured with the same commitment to quality E&J is known for.

Designed for optimal patient transport as well as the hospital/clinical use. Built with a heavy duty frame to support up to 500 lb users. Reinforced at all critical load-bearing points and a long lasting embossed back and seat upholstery with a double inner liner for durability. With its multiple seat widths ranging from 20" to 24" and 18" seat depth (1" longer than standard), this heavy duty wheelchair enables institutions and caregivers to accommodate a large variety of patients.

### Features

- Available in 20", 22" and 24" seat widths, with 18" seat depth (1" longer than standard)
- Carbon steel reinforced frame include double cross braces, heavy-duty 5/8" diameter axles
- Embossed black leatherette upholstery with double inner lining for durability with standard chart pocket on back
- Hammertone paint provides a long lasting and durable finish
- Two-position axle and caster fork provide 2" seat height adjustment for a hemi position
- Tool-free adjustable front riggings with durable aluminum black footplates
- Rear 24" x 1" black plastic wheels with black handrims and non-marking polyurethane tires
- Maximum Weight Capacity 500 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007

### Warranty:

Five (5) year limited warranty



desk arms (shown) and  
swingaway footrest (shown)  
— other configurations available



Anti-fold bar



Wheel-lock



with Elevating Legrests



Color variations may occur due to monitor or printer settings.



## Traveler® HD Heavy Duty Wheelchair

| Item     | HCPCS* code | Overall Chair Width |      | Overall Chair Depth |      | Overall Chair Height |      | Folded Chair Width |      | Seat Width |      | Seat Depth |      | Seat Height |      | Back Height |      | Wheel Diameter |    | Caster Diameter |      | Product <sup>8</sup> Weight w/o front rigging |      | Weight <sup>1</sup> Capacity EVENLY DISTRIBUTED |       | Detach. Full arms | Detach. Desk Fixed Full | Swingaway Footrest | Elevating Legrest | Fixed Footrest | Rear Anti-Tippers |
|----------|-------------|---------------------|------|---------------------|------|----------------------|------|--------------------|------|------------|------|------------|------|-------------|------|-------------|------|----------------|----|-----------------|------|-----------------------------------------------|------|-------------------------------------------------|-------|-------------------|-------------------------|--------------------|-------------------|----------------|-------------------|
|          |             | inch                | cm   | inch                | cm   | inch                 | cm   | inch               | cm   | inch       | cm   | inch       | cm   | inch        | cm   | inch        | cm   | inch           | cm | lb              | kg†  | lb                                            | kg†  |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010320 | K0007       | 28.5                | 72.4 | 34                  | 86.4 | 36                   | 91.4 | 12.5               | 31.8 | 20         | 50.8 | 18         | 45.7 | 18          | 45.7 | 16.5        | 41.9 | 24             | 61 | 8               | 20.3 | 57                                            | 25.9 | 500                                             | 226.8 |                   |                         |                    |                   |                |                   |
| 3G010330 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010340 | K0007       | 28.5                | 72.4 | 34                  | 86.4 | 36                   | 91.4 | 12.5               | 31.8 | 20         | 50.8 | 18         | 45.7 | 20          | 50.8 | 16.5        | 41.9 | 24             | 61 | 8               | 20.3 | 57                                            | 25.9 | 500                                             | 226.8 |                   |                         |                    |                   |                |                   |
| 3G010350 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010420 | K0007       | 30.5                | 77.5 | 34                  | 86.4 | 36                   | 91.4 | 12.5               | 31.8 | 22         | 55.9 | 18         | 45.7 | 20          | 50.8 | 16.5        | 41.9 | 24             | 61 | 8               | 20.3 | 57                                            | 25.9 | 500                                             | 226.8 |                   |                         |                    |                   |                |                   |
| 3G010430 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010440 | K0007       |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010450 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010520 | K0007       | 32.5                | 82.6 | 34                  | 86.4 | 36                   | 91.4 | 12.5               | 31.8 | 24         | 61   | 18         | 45.7 | 20          | 50.8 | 16.5        | 41.9 | 24             | 61 | 8               | 20.3 | 57                                            | 25.9 | 500                                             | 226.8 |                   |                         |                    |                   |                |                   |
| 3G010530 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010540 | K0007       |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G010550 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |
| 3G020500 |             |                     |      |                     |      |                      |      |                    |      |            |      |            |      |             |      |             |      |                |    |                 |      |                                               |      |                                                 |       |                   |                         |                    |                   |                |                   |

Standard

Optional

Available

Accessory

— Not Apply / Available

‡ designed to meet the requirements of

† standard rounding method used

¹ maximum weight capacity - evenly distributed - refer to user manual

⁴ approximate measurement; product weight — as seat width and rigging change, so does product weight

⁵ options available at additional cost - accessories sold separately

■ Standard ■ Optional ■ Available ■ Accessory — Not Apply / Available

‡ designed to meet the requirements of † standard rounding method used † maximum weight capacity - evenly distributed - refer to user manual  
\* approximate measurement: product weight — as seat width and rigging change, so does product weight ★ options available at additional cost - accessories sold separately

### Parts & Accessories: ★ sold separately

|                                                 | HCPCS* | UOM  | Item No. |
|-------------------------------------------------|--------|------|----------|
| Heel loop                                       | E0951  | pair | 90763064 |
| Wheel lock, push-to-lock (includes rubber tips) | E2206  | pair | 90763060 |
| Wheel lock, pull-to-lock (includes rubber tips) | E2206  | pair | 90763061 |
| Wheel lock extension,                           | E0961  | pair | 90763062 |

|                                                      | HCPCS* | UOM     | Item No. |
|------------------------------------------------------|--------|---------|----------|
| Rubber tip, wheel lock                               | —      | 10 each | 90763063 |
| Positioning belt (auto-style) push button 27" to 47" | E0978  | each    | 300-505  |
| Anti-theft locking nuts for front rigging            | —      | 4 each  | 90763630 |
| O2 cylinder holder                                   | E2208  | each    | B10325   |

|                                                       | HCPCS* | UOM  | Item No. |
|-------------------------------------------------------|--------|------|----------|
| Anti-tippers, (includes hardware) for 20" seat height | E0971  | pair | 90763065 |
| Anti-tippers, (includes hardware) for 18" seat height | E0971  | pair | 90763066 |
| Anti-fold bar (chromed) lower frame                   | —      | each | 92523850 |
| Wheelchair Positioning Tray                           | —      | each | EJTRAY18 |



GF Health Products, Inc. is an  
ISO 13485:2016 / MDSAP Certified Company

