

Aluminum Transport Chair

The Everest & Jennings Transport Chair is designed as an affordable patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.



7 stylish color options

Standard Features:

- Available in 17" and 19" seat widths
- Aluminum frame in stylish color options of Red, Blue, Silver, Green, Pink, Aqua and Black
- Fold-down back makes transport and storage easy
- Black padded nylon upholstery is flame retardant
- Carry straps on seat allow for ease of lifting and folding
- Full length padded arms provide user comfort
- Tool-free adjustable detachable swingaway footrests
- Comes standard with 3" webbing black nylon leg strap with velcro adjustable loops to slide over installed footrests
- Black padded nylon leg strap fits 17" and 19" seat width transport chair frames
- Attendant wheel locks on both rear casters for better control and access by companion
- Weighs 20 lb (w/o front rigging)
- Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of Medicare code E1038

Warranty:

One (1) year limited warranty



Aluminum Transport Chair		Dimensions and Weights									
		Seat			Back	Overall				Weight	
		Width	Depth	Height	Height	Depth	Height	Width	Folded Width	Without front rigging	Max Capacity
											
EJ781-1	Lightweight Aluminum - Red	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ782-1	Lightweight Aluminum - Red	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
EJ783-1	Lightweight Aluminum - Silver	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ784-1	Lightweight Aluminum - Silver	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
EJ787-1	Lightweight Aluminum - Blue	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ788-1	Lightweight Aluminum - Blue	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
EJ789-1	Lightweight Aluminum - Green	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ790-1	Lightweight Aluminum - Pink	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ791-1	Lightweight Aluminum - Pink	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
EJ793-1	Lightweight Aluminum - Aqua	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb
EJ785-1	Lightweight Aluminum - Black	17"	16"	19 ¼"	18 ½"	23"	39 ¼"	22 ¼"	10 ¾"	20 lb	250 lb
EJ786-1	Lightweight Aluminum - Black	19"	16"	19 ¼"	18 ½"	23"	39 ¼"	24 ¼"	10 ¾"	20 lb	250 lb

All dimensions are ±1/4"
*Without Front Rigging.

Common Accessories & Parts		HCPCS Code
EJ76-BELT	Positioning belt (plastic buckle) 42" for Transport chairs	E0978