

WHEELCHAIRS





Everest & Jennings®

by graham-field

Value:

Graham-Field is known for a broad product offering at affordable pricing, without sacrifice to quality. Like you, we understand that the cost of a product is more than a price tag. The cost also includes the time and effort required to sell to and service your customers. At Graham-Field we are proud to offer the value found in the Everest & Jennings® line of wheelchairs.

Quality:

The E&J™ name has been synonymous with quality manual folding wheelchairs since 1933. In a time when many manufacturers are producing products at the lowest possible cost, often cutting corners on quality, Graham-Field is standing by its dedication to quality through product design, specification control, testing and inspection. It is easy to see why the E&J Traveler® name is synonymous with quality.

Service:

Our products are often only a day away with centralized distribution located in the metro areas of Atlanta, Georgia; St. Louis, Missouri; and Los Angeles, California. Graham-Field also carries a full line of E&J parts and upholstery to meet your service needs. In addition, our customer service department and technical service teams are always ready to help with your product needs.



Transport Chairs

Aluminum Transport Chair

EJ78X-1 & EJ79X-1 Series

The Everest & Jennings Transport Chair is designed as an affordable short-term patient transfer. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

- 17 and 19 inch (43.18 and 48.26 cm) seat widths [19 inch on Deluxe model only]
- Lightweight aluminum frame available in Aqua, Black, Red, Blue, Silver and Pink [Black and Blue Celeste are available for Deluxe Models only]
- Rear mounted attendant locks on rear casters make it easy for a companion to engage
- Fold-down back makes transport and storage easy
- Black, padded, nylon upholstery is flame retardant
- Tool-free adjustable detachable swingaway footrests
- Comfortable fixed full length arms with padded armrests
- Comes standard with 3 inch (7.62 cm) leg adjustable loops and strap to slide over installed footrests
- Includes carry straps on seat for lifting and folding ease
- 42 inch (106.68 cm) length positioning belt
- Weighs 20 lb (9.07 kg) [without front rigging]
- 250 lb (113.40 kg) weight capacity- EVENLY DISTRIBUTED
- One (1) year limited warranty



Steel Transport Chair

EJ79X-1 Series

- Steel frame is finished in a durable hammertone coating
- Weighs 24 lb (10.89 kg) [without front rigging]
- 250 lb (113.40 kg) weight capacity - EVENLY DISTRIBUTED
- One (1) year limited warranty

E1038*



Deluxe Transport Chair

EJ87X-1 Series

- Handle mounted attendant locks
- 12 inch (55.88 kg) flat-free rear wheels with 8 inch (20.32 kg) swivel front casters
- Weighs 21.38 lb (9.70 kg) [without front rigging]
- 300 lb (136.08 kg) weight capacity - EVENLY DISTRIBUTED
- One (1) year limited warranty



Bariatric Transport

This chair satisfies additional compact transport needs with extra capacity. Handle-mounted attendant wheel locks are easily engaged by a companion.

- Dual crossbrace on gusseted aluminum frame with silver epoxy painted finish
- Available in 24 inch (60.96 cm) width, 18 inch (45.72 cm) seat depth
- 12 inch (30.48 cm) flat-free rear wheels with 8 inch (20.32 cm) swivel front casters
- Removable flip-back desk length padded arms provide user comfort
- Detachable swingaway footrests
- Available with anti-tippers on models after serial number IOS291
- Weighs 31 lb (14.06 kg) [without front rigging]
- 400 lb (181.44 kg) weight capacity - EVENLY DISTRIBUTED
- Three (3) year limited warranty
- Item No. EJ777-3

